



Instant Ice Pack

Weiyou® Instant Ice Pack contains a small pouch of water and some Ammonium Nitrate powder inside the white PVC bag. Weiyou® Instant Ice Pack is ideal for temporary relief of swelling, sprains, strains and muscle aches, Weiyou® Instant Ice Pack can also be used for keeping drink or food cold if needed. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

Product Description

Weiyou® Instant Ice Pack contains a small pouch of water and some Ammonium Nitrate powder inside the white PVC bag. When the patients need external cold therapy, just squeezes the bag until the water pouch inside burst. The whole bag will get cold soon, when the powder and the water mix together, then wrap the ice bag with a cloth and apply to the injury area. Weiyou® Instant Ice Pack is ideal for temporary relief of swelling, sprains, strains and muscle aches, Weiyou® Instant Ice Pack can also be used for keeping drink or food cold if needed.

Specification

Item	Instant Ice Pack
Materials	Ammonium Nitrate, water
Customized Logo	Acceptable
Color	Customized
Sizes	- 12.5*15.5cm,100g
	- 15*20cm, 300g
	- Customized size and weight

Features

- Can be squeezed with one hand
- Be cold immediately when needed



www.nbweiyou.com

Size and weight can be customized
For multiple external medical needs

Ningbo Weiyou Import & Export Co., Ltd.
Phone:+86-13567833512
E-mail:dario@nbweiyou.com









